

LAB

Specimen ID: Acct #: Phone: Rte: Control ID:

## յուլիսիմըմիկինիկինկիսակիմիրիկիկիկիկինուրկ

REFERENCE INTERVAL

Patient Details

DOB:

Age(y/m/d):

Gender:

Patient ID:

Date collected:

Date collected:

Date received:

Date entered:

Date reported:

Physician Details

Ordering:

Referring:

ID:

NPI:

FLAG

UNITS

RESULT

**General Comments & Additional Information** 

Clinical Info: NORMAL

## **Ordered Items**

SYNTHETIC OPIOIDS, SCR w/CONF, U

TESTS

|                            |                                                                                             | A Real Property lives |
|----------------------------|---------------------------------------------------------------------------------------------|-----------------------|
| SYNTHETIC OPIOIDS          |                                                                                             |                       |
| SYNTHETIC OPIOIDS          |                                                                                             |                       |
|                            | Negative                                                                                    | 01                    |
| SYNTHETIC OPIOIDS          |                                                                                             | 22200                 |
|                            | Negative                                                                                    | 01                    |
| KRATOM                     | Negative                                                                                    | 01                    |
| Synthetic op:<br>Screening | ioid analysis includes the following analytes:                                              |                       |
|                            | Analytes                                                                                    |                       |
|                            |                                                                                             |                       |
| SYNTHETIC OP:              | IOIDS GROUP A                                                                               |                       |
| 0 25 ng/mI. 1              | Fentanyl, Acetyl fentanyl                                                                   |                       |
| 0.50 ng/mL 1               |                                                                                             |                       |
|                            | Acetyl norfentanyl                                                                          |                       |
| CAMBITEMIC OD              | IOIDS GROUP B                                                                               |                       |
| SINIHETIC OP.              | TOIDS GROUP B                                                                               |                       |
| 0.25 ng/mL J               | Furanyl fentanyl, p-Fluorobutyrylfentanyl,                                                  |                       |
|                            | Acrylfentanyl, Butyryl fentanyl, Carfentanil,                                               |                       |
|                            | 3-Methyl fentanyl, Valeryl fentanyl,<br>Methoxyacetyl fentanyl, Cyclopropyl fentanyl,       |                       |
|                            | Tetrahydrofuran fentanyl, Phenylacetyl fentanyl,                                            |                       |
|                            | p-Fluoro cyclopropyl fentanyl, MT-45,                                                       |                       |
|                            | Butyryl fentanyl carboxy metabolite,                                                        |                       |
|                            | Isotonitazene, Metonitazene                                                                 |                       |
|                            | Acetylfentanyl-4-methylphenethyl, Ocfentanil, p-Fluorofentanyl, 4-Methoxy-butyryl fentanyl, |                       |
|                            | AH-7921, U-47700, AH-8533, AH-8529, IC-26                                                   |                       |
| 1.00  ng/mL I              | B-Hydroxythiofentanyl, Norcarfentanil,                                                      |                       |
|                            | Xylazine, 4-Hydroxy xylazine                                                                |                       |
| 25.0 ng/mL I               | Methoxyacetyl norfentanyl                                                                   |                       |
| KRATOM                     |                                                                                             |                       |
|                            |                                                                                             |                       |
| 1.00 ng/mL N               | Mitragynine, 7-OH-mitragynine                                                               |                       |

## **Patient Report**



Patient:
DOB: Patient ID: Control ID: Date collected:

TESTS RESULT FLAG UNITS REFERENCE INTERVAL LAB

Analysis performed by Liquid Chromatography with Tandem Mass Spectrometry (LC/MS/MS).

This test was developed and its performance characteristics determined by Labcorp. It has not been cleared or approved by the Food and Drug Administration.