



Patient Information	Specimen Information	Client Information
DOB: AGE:	Specimen: Requisition: Lab Ref #:	
Gender: Phone: Patient ID:	Collected: Received: Reported:	

Test Name In Range Out Of Range Reference Range CALPROTECTIN, STOOL <5 mcg/g EZ

Reference Range: <50 Normal 50-120 Borderline >120 Elevated

Calprotectin in Crohn's disease and ulcerative colitis can be five to several thousand times above the reference population (50 mcg/g or less). Levels are usually 50 mcg/g or less in healthy patients and with irritable bowel syndrome. Repeat testing in 4--6 weeks is suggested for borderline values.