

Specimen ID:
Control ID:

Phone:

Rte:



Patient Details

DOB:
Age(y/m/d):
Gender: SSN:
Patient ID:

Specimen Details

Date collected:
Date received:
Date entered:
Date reported:

Physician Details

Ordering:
Referring:
ID:
NPI:

General Comments & Additional Information

Alternate Control Number:
Total Volume:

Alternate Patient ID:
Fasting:

Ordered Items

Metanephrines, Frac., Pl. Free

TESTS	RESULT	FLAG	UNITS	REFERENCE INTERVAL	LAB
Metanephrines, Frac., Pl. Free					
Normetanephrine, Pl	86		pg/mL	0 - 145	01
Metanephrine, Pl	37		pg/mL	0 - 62	01

Concentrations of Normetanephrine between 146 and 487 pg/mL, and Metanephrine between 63 and 255 pg/mL are considered indeterminate. Follow-up biochemical testing is recommended when patient levels fall within this indeterminate range. These tests include repeat testing of plasma/urinary fractionated metanephrines and plasma catecholamines.