

Specimen ID:
 Control ID:

Acct #:

Phone:

Rte:


Patient Details

 DOB:
 Age(y/m/d):
 Gender: SSN:
 Patient ID:

Specimen Details

 Date collected:
 Date received:
 Date entered:
 Date reported:

Physician Details

 Ordering:
 Referring:
 ID:
 NPI:

General Comments & Additional Information

Clinical Info: NORMAL REPORT

Total Volume: Not Provided

Fasting: No

Ordered Items

Nickel, Blood

TESTS	RESULT	FLAG	UNITS	REFERENCE INTERVAL	LAB
Nickel, Blood					
Nickel	<2.0		ng/mL	<2.0	01
Serum is the preferred specimen for monitoring nickel exposure. Nickel analysis performed by inductively coupled plasma / mass spectrometry (ICP/MS).					

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.